



Product Specification
S-Monovette® Serum CAT, 5.5 ml, cap white,
(LxØ): 75 x 15 mm, with paper label

Page 1



Product description

| | |
|---------------------|--|
| Order number | 03.1397.001 |
| Product description | S-Monovette® Serum CAT, preparation: Clotting activator, 5.5 ml, membrane screw cap, cap white, colour code EU, (LxØ) without cap: 75 x 15 mm, with paper label, label/print: black, 50 piece(s)/case, sterile |

Product characteristics

| | |
|----------------------------|--------------------|
| Type of collection | venous |
| Type of preparation | Clotting activator |
| Preparation concentration | 1.6 mg/ml blood |
| Colour code | EU |
| Label/ Print | with paper label |
| Colour of print/label | black |
| Graduation | fill mark |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Recommended centrifugation | 2000 x g |
| Centrifugation temperature | 20 |
| Centrifugation time | 10 Min |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® Serum CAT, 5.5 ml, cap white,
(LxØ): 75 x 15 mm, with paper label

Page 3

| | |
|-------------------|---------------|
| Width of box | 117 mm |
| Height of box | 82 mm |
| Depth of case | 440 mm |
| Width of case | 400 mm |
| Height of case | 129 mm |
| Case volume | 0.0227 cbm |
| Weight of product | 0.007 kg |
| Weight of case | 4.28 kg |
| EAN of inner box | 4038917128438 |
| EAN case | 4038917021722 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT