



Product description

| | |
|---------------------|--|
| Order number | 03.1628 |
| Product description | S-Monovette® Lithium heparin LH, preparation: Lithium heparin, 5.5 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 75 x 15 mm, with plastic label, label/print: transparent/orange, 50 piece(s)/case, sterile |

Product characteristics

| | |
|----------------------------|--------------------|
| Type of collection | venous |
| Type of preparation | Lithium heparin |
| Preparation concentration | 16 IU/ml blood |
| Colour code | EU |
| Label/ Print | with plastic label |
| Colour of print/label | transparent/orange |
| Graduation | fill mark |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Recommended centrifugation | 2000 x g |
| Centrifugation temperature | 20 |
| Centrifugation time | 10 Min |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin LH, 5.5 ml, cap orange, (LxØ): 75 x 15 mm, with plastic label

Page 2

Liquid-tight yes
Connection membrane screw cap

Size

Sample volume 5.5 ml
Outer diameter of cap 18 mm
Diameter 15 mm
Length including cap 91 mm
Length excluding cap 75 mm
Length including cap and plunger 112 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap orange
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest case
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 34000
Depth of box 189 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin LH, 5.5 ml, cap orange, (LxØ): 75 x 15 mm, with plastic label

| | |
|-------------------|---------------|
| Width of box | 117 mm |
| Height of box | 82 mm |
| Depth of case | 440 mm |
| Width of case | 400 mm |
| Height of case | 129 mm |
| Case volume | 0.0227 cbm |
| Weight of product | 0.0069 kg |
| Weight of case | 4.27 kg |
| EAN of inner box | 4038917087285 |
| EAN case | 4038917005746 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

