



Product description

Order number	03.1631.001
Product description	S-Monovette® Lithium heparin gel LH, preparation: Lithium heparin gel, 4.7 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 75 x 15 mm, with paper label, label/print: orange, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin gel
Preparation concentration	25 IU/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	orange
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	15 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin gel LH, 4.7 ml, cap orange, (LxØ): 75 x 15 mm, with paper label

Page 2

Liquid-tight yes
Connection membrane screw cap

Size

Sample volume 4.7 ml
Outer diameter of cap 18 mm
Diameter 15 mm
Length including cap 91 mm
Length excluding cap 75 mm
Length including cap and plunger 112 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap orange
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest case
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 26000
Depth of box 189 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin gel LH, 4.7 ml, cap orange, (LxØ): 75 x 15 mm, with paper label

Page 3

Width of box	117 mm
Height of box	82 mm
Depth of case	440 mm
Width of case	400 mm
Height of case	129 mm
Case volume	0.0227 cbm
Weight of product	0.0075 kg
Weight of case	4.54 kg
EAN of inner box	4038917132206
EAN case	4038917024310

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT