



Product Specification
S-Monovette® Serum CAT, 2.6 ml, cap white,
(LxØ): 65 x 13 mm, with paper label

Page 1



Product description

| | |
|---------------------|--|
| Order number | 04.1904.001 |
| Product description | S-Monovette® Serum CAT, preparation: Clotting activator, 2.6 ml, membrane screw cap, cap white, colour code EU, (LxØ) without cap: 65 x 13 mm, with paper label, label/print: black, 50 piece(s)/case, sterile |

Product characteristics

| | |
|----------------------------|--------------------|
| Type of collection | venous |
| Type of preparation | Clotting activator |
| Colour code | EU |
| Label/ Print | with paper label |
| Colour of print/label | black |
| Graduation | fill mark |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Recommended centrifugation | 2000 x g |
| Centrifugation temperature | 20 |
| Centrifugation time | 10 Min |
| Liquid-tight | yes |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® Serum CAT, 2.6 ml, cap white,
(LxØ): 65 x 13 mm, with paper label

Connection membrane screw cap

Size

Sample volume 2.6 ml
Outer diameter of cap 16 mm
Diameter 13 mm
Length including cap 81 mm
Length excluding cap 65 mm
Length including cap and plunger 102 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap white
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest case
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 42000
Depth of box 171 mm
Width of box 105 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Serum CAT, 2.6 ml, cap white,
(LxØ): 65 x 13 mm, with paper label

| | |
|-------------------|---------------|
| Height of box | 73 mm |
| Depth of case | 388 mm |
| Width of case | 358 mm |
| Height of case | 125 mm |
| Case volume | 0.0174 cbm |
| Weight of product | 0.0043 kg |
| Weight of case | 2.87 kg |
| EAN of inner box | 4038917134224 |
| EAN case | 4038917025539 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

