



Product Specification

S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 3 ml, cap blue, (LxØ): 75 x 13 mm, with paper label

Page 1



Product description

| | |
|---------------------|---|
| Order number | 04.1919.100 |
| Product description | S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, preparation: Citrate 9NC, 3 ml, membrane screw cap, cap blue, colour code ISO, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: blue, 50 piece(s)/case, sterile |

Product characteristics

| | |
|----------------------------|--------------------|
| Type of collection | venous |
| Type of preparation | Citrate 9NC |
| Preparation concentration | 0.106 mol/l |
| Colour code | ISO |
| Label/ Print | with paper label |
| Colour of print/label | blue |
| Graduation | fill mark |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Recommended centrifugation | 1800 x g |
| Centrifugation temperature | 20 |
| Centrifugation time | 10 Min |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 3 ml,
cap blue, (LxØ): 75 x 13 mm, with paper label**

Page 3

| | |
|-------------------|---------------|
| Width of box | 116 mm |
| Height of box | 73 mm |
| Depth of case | 388 mm |
| Width of case | 358 mm |
| Height of case | 135 mm |
| Case volume | 0.0188 cbm |
| Weight of product | 0.0045 kg |
| Weight of case | 3.01 kg |
| EAN of inner box | 4038917266970 |
| EAN case | 4038917266963 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT