



## Product Specification

**S-Monovette® Lithium heparin LH, fluid, 4.9 ml, cap orange, (LxØ): 90 x 13 mm, with paper label**

**Page 1**



### Product description

Order number	04.1939.001
Product description	S-Monovette® Lithium heparin LH, fluid, preparation: Lithium-Heparin, liquid, 4.9 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 90 x 13 mm, with paper label, label/print: orange, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Lithium-Heparin, liquid
Preparation concentration	19 IU/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	orange
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**





## Product Specification

**S-Monovette® Lithium heparin LH, fluid, 4.9 ml,  
cap orange, (LxØ): 90 x 13 mm, with paper label**

**Page 3**

Width of box	130 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	150 mm
Case volume	0.0208 cbm
Weight of product	0.0048 kg
Weight of case	3.18 kg
EAN of inner box	4038917134781
EAN case	4038917053761

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**