



## Product Specification

**S-Monovette® GlucoEXACT FC, 3.1 ml, cap pink,  
(LxØ): 75 x 13 mm, with paper label**

**Page 1**

### Product description

Order number	04.1945.201
Product description	S-Monovette® GlucoEXACT FC, preparation: Citrate/fluoride/EDTA, 3.1 ml, membrane screw cap, cap pink, colour code ISO, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: pink, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Citrate/fluoride/EDTA
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	pink
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes
Connection	membrane screw cap

### Size

Sample volume	3.1 ml
Outer diameter of cap	16 mm
Diameter	13 mm
Length including cap	91 mm
Length excluding cap	75 mm
Length including cap and plunger	112 mm

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	pink
Membrane material	rubber
Tube material	Polypropylene (PP)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® GlucoEXACT FC, 3.1 ml, cap pink,  
(LxØ): 75 x 13 mm, with paper label**

**Page 2**

Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

### Purity & certification

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

### Packaging

Minimum order qty. 20000  
Type of smallest case  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 39000  
Depth of box 171 mm  
Width of box 116 mm  
Height of box 73 mm  
Depth of case 388 mm  
Width of case 358 mm  
Height of case 135 mm  
Case volume 0.0188 cbm  
Weight of product 0.0047 kg  
Weight of case 3.08 kg  
EAN of inner box 4038917554091  
EAN case 4038917554084

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**