



Product Specification

**S-Monovette® Lithium heparin gel+ LH, 2.7 ml,
cap orange, (LxØ): 75 x 13 mm, with paper label**

Page 1



Product description

Order number	04.1952.001
Product description	S-Monovette® Lithium heparin gel+ LH, preparation: Lithium heparin gel+, 2.7 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: brown, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin gel+
Preparation concentration	25 IU/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	brown
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	7 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**S-Monovette® Lithium heparin gel+ LH, 2.7 ml,
cap orange, (LxØ): 75 x 13 mm, with paper label**

Page 3

Width of box	116 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0051 kg
Weight of case	3.31 kg
EAN of inner box	4038917523035
EAN case	4038917523028

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT