



## Product Specification

**S-Monovette® Lithium heparin gel+ LH, 2.7 ml, cap orange, (LxØ): 75 x 13 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix**

Page 1



### Product description

Order number	04.1952.088
Product description	S-Monovette® Lithium heparin gel+ LH, preparation: Lithium heparin gel+, 2.7 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 75 x 13 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix 034, label/print: white/brown, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin gel+
Preparation concentration	25 IU/ml blood
Colour code	EU
Label/ Print	with plastic label pre-barcoded
Colour of print/label	white/brown
Graduation	fill mark
Cap	membrane screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

Page 2

**S-Monovette® Lithium heparin gel+ LH, 2.7 ml, cap orange, (LxØ): 75 x 13 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix**

Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	7 Min
Liquid-tight	yes
Connection	membrane screw cap

### Size

Sample volume	2.7 ml
Outer diameter of cap	16 mm
Diameter	13 mm
Length including cap	91 mm
Length excluding cap	75 mm
Length including cap and plunger	112 mm

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	orange
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Lithium heparin gel+ LH, 2.7 ml, cap orange, (LxØ): 75 x 13 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix**

**Page 3**

### Packaging

Minimum order qty.	20000
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	39000
Depth of box	171 mm
Width of box	116 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0051 kg
Weight of case	3.29 kg
EAN of inner box	4038917950060
EAN case	4038917950077

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**