



## Product Specification

**S-Monovette® Lithium heparin gel+ LH, 4.9 ml, cap orange, (LxØ): 90 x 13 mm, with plastic label**

**Page 1**



### Product description

Order number	04.1954
Product description	S-Monovette® Lithium heparin gel+ LH, preparation: Lithium heparin gel+, 4.9 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 90 x 13 mm, with plastic label, label/print: transparent/brown, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin gel+
Preparation concentration	25 IU/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/brown
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Lithium heparin gel+ LH, 4.9 ml,  
cap orange, (LxØ): 90 x 13 mm, with plastic label**

Page 2

Centrifugation time	7 Min
Liquid-tight	yes
Connection	membrane screw cap

### Size

Sample volume	4.9 ml
Outer diameter of cap	16 mm
Diameter	13 mm
Length including cap	106 mm
Length excluding cap	90 mm
Length including cap and plunger	127 mm

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	orange
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	36000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Lithium heparin gel+ LH, 4.9 ml,  
cap orange, (LxØ): 90 x 13 mm, with plastic label**

**Page 3**

Depth of box	171 mm
Width of box	130 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	150 mm
Case volume	0.0208 cbm
Weight of product	0.0057 kg
Weight of case	3.61 kg
EAN of inner box	4038917420327
EAN case	4038917420310

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**