



## Product description

Order number	05.1073.100
Product description	S-Monovette® Fluoride/EDTA FE, preparation: Fluoride + EDTA, 2.7 ml, membrane screw cap, cap grey, colour code ISO, (LxØ) without cap: 66 x 11 mm, with paper label, label/print: grey, 50 piece(s)/case, sterile

## Product characteristics

Type of collection	venous
Type of preparation	Fluoride + EDTA
Preparation concentration	1.0 mg/ml blood + 1.2 mg/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	grey
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	9 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.







## Product Specification

**S-Monovette® Fluoride/EDTA FE, 2.7 ml, cap grey,  
(LxØ): 66 x 11 mm, with paper label**

**Page 3**

Width of box	108 mm
Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0039 kg
Weight of case	2.48 kg
EAN of inner box	4038917125451
EAN case	4038917020060

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**