



### Product description

Order number	05.1165.100
Product description	S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, preparation: Citrate 9NC, 3 ml, membrane screw cap, cap blue, colour code ISO, (LxØ) without cap: 66 x 11 mm, with paper label, label/print: blue, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Citrate 9NC
Preparation concentration	0.106 mol/l
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	blue
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	1800 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 3 ml,  
cap blue, (LxØ): 66 x 11 mm, with paper label**

**Page 2**

Liquid-tight Connection yes  
membrane screw cap

### Size

Sample volume 3 ml  
Outer diameter of cap 13 mm  
Diameter 11 mm  
Length including cap 82 mm  
Length excluding cap 66 mm  
Length including cap and plunger 103 mm

### Material & colours

Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap blue  
Membrane material rubber  
Tube material Polypropylene (PP)  
Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

### Purity & certification

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

### Packaging

Minimum order qty. 500  
Type of smallest subpackaging case  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 56000  
Depth of box 145 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 3 ml,  
cap blue, (LxØ): 66 x 11 mm, with paper label**

**Page 3**

Width of box	108 mm
Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0041 kg
Weight of case	2.61 kg
EAN of inner box	4038917126113
EAN case	4038917020442

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**