




Product Specification
**S-Monovette® neutral Z, 2.7 ml, cap white, (LxØ):
66 x 11 mm, with plastic label**

Page 1



Product description

Order number	05.1729
Product description	S-Monovette® neutral Z, without preparation, 2.7 ml, membrane screw cap, cap white, colour code ISO, (LxØ) without cap: 66 x 11 mm, with plastic label, label/print: white, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	without preparation
Colour code	ISO
Label/ Print	with plastic label
Colour of print/label	white
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Coding colour	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® neutral Z, 2.7 ml, cap white, (LxØ):
66 x 11 mm, with plastic label

Page 3

Width of box	108 mm
Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0038 kg
Weight of case	2.43 kg
EAN of inner box	4038917240895
EAN case	4038917240888

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT