



**Product Specification**  
**S-Monovette® EDTA K2E, 1.2 ml, cap violet,**  
**(LxØ): 66 x 8 mm, with plastic label**

**Page 1**



### Product description

Order number	06.1660.100
Product description	S-Monovette® EDTA K2E, preparation: K2 EDTA, 1.2 ml, membrane screw cap, cap violet, colour code ISO, (LxØ) without cap: 66 x 8 mm, with plastic label, label/print: white/violet, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	K2 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	ISO
Label/ Print	with plastic label
Colour of print/label	white/violet
Graduation	yes
Cap	membrane screw cap
Closure type	screw cap
Liquid-tight	yes
Connection	membrane screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® EDTA K2E, 1.2 ml, cap violet,**  
**(LxØ): 66 x 8 mm, with plastic label**

**Page 2**

**Size**

Sample volume	1.2 ml
Outer diameter of cap	13 mm
Diameter	8 mm
Length including cap	82 mm
Length excluding cap	66 mm
Length including cap and plunger	103 mm

**Material & colours**

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	violet
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	rubber
Piston plunger material	Polystyrene (PS)

**Purity & certification**

Satisfies the requirement	IATA, ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	56000
Depth of box	145 mm
Width of box	108 mm
Height of box	63 mm
Depth of case	333 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® EDTA K2E, 1.2 ml, cap violet,**  
**(LxØ): 66 x 8 mm, with plastic label**

**Page 3**

Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0033 kg
Weight of case	2.18 kg
EAN of inner box	4038917322737
EAN case	4038917322720

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**