



Product description

Order number	65.386.992
Product description	Colour-coded inserts, for CryoPure tubes, colour mix, 100 piece(s)/bag

Product characteristics

Graduation	no
Suitable for	CryoPure tubes

Material & colours

Product material	Polystyrene (PS)
Colour of product	colour mix

Purity & certification

Product category	no medical device IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	2500
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	2500
Piece(s) / pallet	810000
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.0002 kg
Weight of case	0.6 kg
EAN of inner box	4038917228541
EAN case	4038917228534

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

