




**Product Specification**  
**Tissue culture dish, (ØxH): 35 x 10 mm, surface:**  
**Standard**

**Page 1**



### Product description

Order number	83.3900.002
Product description	Tissue culture dish, (ØxH): 35 x 10 mm, with grid, material: PS, surface: Standard, for Adherent cells, coding colour: red, TC Tested, 10 piece(s)/bag

### Product characteristics

Closure type	slip-on lid
Cell type	Adherent cells
Base shape	SUREGrip Design
Growth area	9.4 cm <sup>2</sup>
Volume of work	3 ml
Vessel type	Dish
Coding colour	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Tissue culture dish, (ØxH): 35 x 10 mm, surface:**  
**Standard**

**Page 2**

**Size**

Diameter	35 mm
Height of product	10 mm

**Material & colours**

Product material	Polystyrene (PS)
Closure material	Polystyrene (PS)

**Purity & certification**

Product category	no medical device   IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	10
Piece (s) / outer case	500
Piece(s) / pallet	30000
Depth of case	322 mm
Width of case	257 mm
Height of case	245 mm
Case volume	0.0203 cbm
Weight of product	0.0042 kg
Weight of case	2.82 kg
EAN of inner box	4038917337984
EAN case	4038917365796

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**