



Product description

Order number	83.3900.300
Product description	Tissue culture dish, (ØxH): 35 x 10 mm, material: PS, surface: Cell+, for challenging adherent cells, coding colour: yellow, TC Tested, 10 piece(s)/bag

Product characteristics

Closure type	slip-on lid
Cell type	challenging adherent cells
Base shape	SUREGrip Design
Growth area	9.4 cm ²
Volume of work	3 ml
Vessel type	Dish
Coding colour	yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tissue culture dish, (ØxH): 35 x 10 mm, surface:
Cell+

Page 2

Size

Diameter	35 mm
Height of product	10 mm

Material & colours

Product material	Polystyrene (PS)
Closure material	Polystyrene (PS)

Purity & certification

Product category	no medical device IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	10
Piece (s) / outer case	500
Piece(s) / pallet	30000
Depth of case	322 mm
Width of case	257 mm
Height of case	245 mm
Case volume	0.0203 cbm
Weight of product	0.0037 kg
Weight of case	2.56 kg
EAN of inner box	4038917338042
EAN case	4038917365802

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT