



Product description

Order number	83.3901.002
Product description	Tissue culture dish, (ØxH): 60 x 15 mm, with grid, material: PS, surface: Standard, for Adherent cells, coding colour: red, TC Tested, 10 piece(s)/bag

Product characteristics

Closure type	slip-on lid
Cell type	Adherent cells
Base shape	SUREGrip Design
Growth area	22.1 cm ²
Volume of work	5 ml
Vessel type	Dish
Coding colour	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tissue culture dish, (ØxH): 60 x 15 mm, surface:
Standard

Page 2

Size

Diameter	60 mm
Height of product	15 mm

Material & colours

Product material	Polystyrene (PS)
Closure material	Polystyrene (PS)

Purity & certification

Product category	no medical device IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	10
Piece (s) / outer case	500
Piece(s) / pallet	12000
Depth of case	332 mm
Width of case	332 mm
Height of case	380 mm
Case volume	0.0419 cbm
Weight of product	0.0059 kg
Weight of case	4.02 kg
EAN of inner box	4038917338103
EAN case	4038917338097

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-08-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT