



### **Product description**

Order number	83.3910.300
Product description	Cell culture flask, T-25, material: PS, surface: Cell+, for challenging adherent cells, coding colour: yellow, 2-position screw cap, TC Tested, 10 piece(s)/bag

### **Product characteristics**

Cap	2-position screw cap
Cell type	challenging adherent cells
Growth area	25 cm <sup>2</sup>
Volume of work	7 ml - 12,5 ml
Vessel type	flask
Coding colour	yellow

### **Material & colours**

Product material	Polystyrene (PS)
Closure material	Polyethylene (PE)
Colour of cap	yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category	no medical device   IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

### **Packaging**

Minimum order qty.	300
Type of smallest subpackaging	bag
Piece(s) / inner box	10
Piece (s) / outer case	300
Piece(s) / pallet	6000
Depth of case	602 mm
Width of case	408 mm
Height of case	340 mm
Case volume	0.0835 cbm
Weight of product	0.0152 kg
Weight of case	7.42 kg
EAN of inner box	4038917337410
EAN case	4038917337403

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

