



### Product description

Order number	83.3930.041
Product description	TC insert, for 6-well plates, membrane: PET, transparent, pore size: 0.4 $\mu\text{m}$ , sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

### Product characteristics

Suitable for	6-well plates
Base shape	flat base
Growth area	4.5 cm <sup>2</sup>
Volume of work	1 ml - 4 ml
Vessel type	TC insert
Pore size	0,4

### Size

Height of product	16.38 mm
-------------------	----------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Membrane material	Polyethylene terephthalate (PET)
Frame material	Polystyrene (PS)

### **Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	X-Ray
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

### **Packaging**

Minimum order qty.	24
Type of smallest subpackaging	blister
Piece(s) / inner box	24
Piece (s) / outer case	720
Piece(s) / pallet	11520
Depth of box	196 mm
Width of box	104 mm
Height of box	131 mm
Depth of case	410 mm
Width of case	435 mm
Height of case	550 mm
Case volume	0.0981 cbm
Weight of product	0.0051 kg
Weight of case	6.46 kg
EAN of inner box	4038917370868
EAN case	4038917370851

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

