



Product description

Order number	83.3941.505
Product description	Collection vessel, for Filtropur BT 50, 500 ml, screw cap, cap assembled, sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase, 1 piece(s)/bag

Product characteristics

Graduation	yes
Closure type	screw cap
Suitable for	Filtropur BT 50
Volume of work	500 ml
Vessel type	flask
Filter	no

Size

Sample volume	500 ml
Filtration volume	500 ml
Height of product	13.93 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase
Sterilisation	Electron irradiation
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	24
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	24
Piece(s) / pallet	576
Depth of case	512 mm
Width of case	391 mm
Height of case	329 mm
Case volume	0.0659 cbm
Weight of product	0.0096 kg
Weight of case	2.14 kg
EAN of inner box	4038917968911
EAN case	4038917968928

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

