



Product description

Order number	83.3941.511
Product description	Filtropur BT 50, Bottle top filter, for screw cap 45 mm, 500 ml, pore size: 0.22 µm, for sterile filtration, sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase, 1 piece(s)/bag

Product characteristics

Graduation	yes
Example of use	for sterile filtration
Suitable for	screw cap 45 mm
Filter	yes
Surface area of filter	38.04 cm ²
Pore size	0,22
Connection	Step taper

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Sample volume	500 ml
Filtration volume	500 ml
Height of product	12.87 cm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Filter material	Polyethersulfone (PES)

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase
Sterilisation	Electron irradiation
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	24
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	24
Piece(s) / pallet	576
Depth of case	510 mm
Width of case	391 mm
Height of case	329 mm
Case volume	0.0668 cbm
Weight of product	0.0684 kg
Weight of case	3.56 kg
EAN of inner box	4038917969024
EAN case	4038917969031

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

