



Product description

Order number	83.3944
Product description	Pressure filtration unit, Filtropur L 0.2 S, pre-filter: GF, membrane: PES, filtration surface: 20 cm ² , pore size: 0.2 µm, for sterile filtration, connection: Tube connector, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

Product characteristics

Example of use	for sterile filtration
Filter	yes
Surface area of filter	20 cm ²
Pore size	0,2
Connection	Tube connector

Size

Filtration volume	10000 ml
-------------------	----------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polycarbonate (PC)
Membrane material	Polyethersulfone (PES)

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gas sterilisation
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	10
Type of smallest subpackaging	bag
Piece(s) / inner box	10
Piece (s) / outer case	50
Piece(s) / pallet	400
Depth of box	352 mm
Width of box	155 mm
Height of box	123 mm
Depth of case	530 mm
Width of case	370 mm
Height of case	280 mm
Case volume	0.0549 cbm
Weight of product	0.025 kg
Weight of case	4.24 kg
EAN of inner box	4038917926225
EAN case	4038917926232

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

