



Product description

Order number 94.6077.037
Product description miniPERM® , Screw cap, for miniPERM® nutrient module, not autoclavable, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Product characteristics

Suitable for miniPERM® nutrient module

Material & colours

Product material High Density Polyethylene (HD-PE)

Purity & certification

Product category no medical device | IVD
Purity standard sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation Gamma irradiation
Pyrogen-free yes
Non-cytotoxic yes
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	16
Type of smallest subpackaging	blister
Piece(s) / inner box	16
Piece (s) / outer case	64
Piece(s) / pallet	2304
Depth of case	388 mm
Width of case	253 mm
Height of case	225 mm
Case volume	0.0221 cbm
Weight of product	0.0081 kg
Weight of case	1.2 kg
EAN of inner box	4038917255776
EAN case	4038917255769

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

