



### **Product description**

Order number	94.6077.135
Product description	miniPERM <sup>®</sup> , Luer cannula, for miniPERM <sup>®</sup> bioreactor, Luer (f), not autoclavable, sterile, non-pyrogenic/endotoxin-free, 1 piece(s)/blister

### **Product characteristics**

Suitable for	miniPERM <sup>®</sup> bioreactor
Connection	Luer (f)

### **Material & colours**

Product material	stainless steel
------------------	-----------------

### **Purity & certification**

Product category	Medical device
CE certified	CE supplier
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	100
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	1200
Piece(s) / pallet	129600
Depth of case	383 mm
Width of case	193 mm
Height of case	190 mm
Case volume	0.014 cbm
Weight of product	0.0011 kg
Weight of case	1.8 kg
EAN of inner box	4038917255899
EAN case	4038917255882

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

