



**Product Specification**  
**miniPERM<sup>®</sup> , Disposable syringe 50 ml, for  
miniPERM<sup>®</sup> bioreactor, Luer lock**

**Page 1**



**Product description**

Order number	94.6077.137
Product description	miniPERM <sup>®</sup> , Disposable syringe 50 ml, for miniPERM <sup>®</sup> bioreactor, Luer lock, not autoclavable, sterile, 1 piece(s)/blister

**Product characteristics**

Suitable for	miniPERM <sup>®</sup> bioreactor
Connection	Luer lock

**Size**

Max. volume	50 ml
-------------	-------

**Material & colours**

Product material	Polypropylene (PP)
------------------	--------------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**miniPERM® , Disposable syringe 50 ml, for**  
**miniPERM® bioreactor, Luer lock**

**Page 2**

**Purity & certification**

Product category	In vitro diagnostic device, CE
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	60
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	60
Piece(s) / pallet	3840
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0.0255 cbm
Weight of product	0.0366 kg
Weight of case	2.53 kg
EAN of inner box	4038917255950
EAN case	4038917255943

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**