



Product description

Order number	94.6077.138
Product description	miniPERM® , Filling tube 5", for miniPERM® bioreactor, not autoclavable, sterile, 1 piece(s)/blister

Product characteristics

Suitable for	miniPERM® bioreactor
--------------	----------------------

Size

Diameter	5.9 mm
Width of product	5.9 mm
Height of product	127 mm
Length of product	5.9 mm

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**miniPERM® , Filling tube 5" , for miniPERM®
bioreactor**

Page 2

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	50
Type of smallest subpackaging	blister
Piece(s) / inner box	50
Piece (s) / outer case	300
Piece(s) / pallet	14400
Depth of case	353 mm
Width of case	268 mm
Height of case	260 mm
Case volume	0.0246 cbm
Weight of product	0.001 kg
Weight of case	1.55 kg
EAN of inner box	4038917376488
EAN case	4038917255974

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT