



### **Product description**

Order number	94.6140.002
Product description	x-well cell culture chamber, flask, on PCA slide, removable frame, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 6 piece(s)/blister

### **Product characteristics**

Closure type	screw cap
Design	on PCA slide
Cell type	Adherent cells
Number of wells	1
Growth area	9 cm <sup>2</sup>
Volume of work	4 ml
Removable frame	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Width of product	24 mm
Length of product	76 mm

### **Material & colours**

Product material	Polycycloalkanes (PCA)
Frame material	Polystyrene (PS)

### **Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

### **Packaging**

Minimum order qty.	6
Type of smallest subpackaging	blister
Piece(s) / inner box	6
Piece (s) / outer case	96
Piece(s) / pallet	8064
Depth of box	127 mm
Width of box	28 mm
Height of box	259 mm
Depth of case	275 mm
Width of case	270 mm
Height of case	255 mm
Case volume	0.0189 cbm
Weight of product	0.0152 kg
Weight of case	2.44 kg
EAN of inner box	4038917351935
EAN case	4038917351928

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

