



### **Product description**

Order number 90.183.100  
Product description Tube holder, SC 2700, for centrifuge SC 2700

### **Product characteristics**

Suitable for centrifuge SC 2700

### **Size**

Diameter 18.4 mm  
Width of product 23 mm  
Height of product 83.5 mm  
Length of product 23 mm

### **Material & colours**

Colour of product black

### **Purity & certification**

Product category no medical device | IVD  
Batched no

### **Packaging**

Minimum order qty. 2  
Type of smallest subpackaging bag  
Piece(s) / inner box 2  
Piece (s) / outer case 50  
Piece(s) / pallet 16200  
Depth of case 268 mm  
Width of case 191 mm  
Height of case 100 mm  
Case volume 0.0051 cbm  
Weight of product 0.021 kg  
Weight of case 1.2 kg  
EAN of inner box 4038917909143  
EAN case 4038917909150

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

