



Product description

Order number 90.183.200
Product description adapter, for centrifuge SC 2700

Product characteristics

Suitable for centrifuge SC 2700

Size

Diameter 17 mm
Width of product 17 mm
Height of product 12.5 mm
Length of product 17 mm

Material & colours

Product material Polypropylene (PP)
Colour of product grey

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 2
Type of smallest subpackaging bag
Piece(s) / inner box 2
Piece (s) / outer case 300
Piece(s) / pallet 97200
Depth of case 268 mm
Width of case 191 mm
Height of case 100 mm
Case volume 0.0051 cbm
Weight of product 0.002 kg
Weight of case 1 kg
EAN of inner box 4038917909679
EAN case 4038917909686

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

