



**Product Specification**  
**Centrifuge, SMC 6 V, supply voltage: 115 V,  
swing-out rotor 6-fold**

**Page 1**



### Product description

Order number	90.184.401
Product description	Centrifuge, SMC 6 V, supply voltage: 115 V, incl. swing-out rotor 6-fold and tube holder 1-fold for tubes up to 100 mm length, for separation of serum/plasma from the blood cells, swing-out rotor 6-fold

### Product characteristics

Design	swing-out rotor 6-fold
Suitable for	separation of serum/plasma from the blood cells
Number of positions	6
Supply voltage	115 V

### Purity & certification

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## **Packaging**

Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	36
Weight of product	5.7995 kg
Weight of case	5.8 kg
EAN of inner box	4038917534420
EAN case	4038917534413

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

