




Product Specification

Multi-purpose container, max. volume: 500 ml,
(LxØ): 150 x 70 mm, graduated, PP, translucent,
with paper label

Page 1



Product description

Order number	51.9922.815
Product description	Multi-purpose container, max. volume: 500 ml, (LxØ): 150 x 70 mm, translucent, cap: natural, graduated, material: PP, screw cap, cap assembled, with paper label, label/print: black, sterile, 80 piece(s)/bag

Product characteristics

Label/ Print	with paper label
Colour of print/label	black
Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Application	Container for water samples
Base shape	flat base
Volume of work	500 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Multi-purpose container, max. volume: 500 ml,
(LxØ): 150 x 70 mm, graduated, PP, translucent,
with paper label**

Page 2

Size

Max. volume	500 ml
Diameter	70 mm
Height of product	150 mm
Length of product	150 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	translucent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural
Tube material	Polypropylene (PP)

Purity & certification

Satisfies the requirement	sealed at 9"hg/10min
Latex-free	yes
Purity standard	sterile
Pyrogen-free	no
Safety data sheet	no

Packaging

Minimum order qty.	105
Type of smallest subpackaging	bag
Piece(s) / inner box	105
Piece (s) / outer case	105
Piece(s) / pallet	2520
Depth of case	534 mm
Width of case	384 mm
Height of case	519 mm
Case volume	0.1025 cbm
Weight of case	7.35 kg
EAN of inner box	4038917153348
EAN case	4038917234863

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT