



### **Product description**

|                     |   |
|---------------------|---|
| Order number        | 80.630  |
| Product description | Forensic swab, round, in the tube with ventilation membrane, with safety label, ETO treated, ISO 18385, passive drying, stem length: 85 mm, stem material: PS, swab tip round, swab material: viscose, tube: PP, 50 piece(s)/case |

### **Product characteristics**

|              |                                       |
|--------------|---------------------------------------|
| Label/ Print | with safety label                     |
| Cap          | standard cap                          |
| Closure type | slip-on lid                           |
| Design       | in the tube with ventilation membrane |
| Application  | Removal of forensic samples           |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### Size

|                      |        |
|----------------------|--------|
| Length of product    | 175 mm |
| Length excluding cap | 105 mm |
| Length of stem       | 85 mm  |

### Material & colours

|                   |                    |
|-------------------|--------------------|
| Colour of product | white              |
| Closure material  | Polypropylene (PP) |
| Colour of cap     | white              |
| Swab material     | viscose            |
| Stem material     | Polystyrene (PS)   |
| Tube material     | Polypropylene (PP) |

### Purity & certification

|                           |                         |
|---------------------------|-------------------------|
| Satisfies the requirement | ISO 18385               |
| Product category          | no medical device   IVD |
| Purity standard           | ETO treated             |
| Sterilisation             | Gas sterilisation       |
| Batched                   | yes                     |

### Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 500           |
| Type of smallest subpackaging | case          |
| Piece(s) / inner box          | 50            |
| Piece (s) / outer case        | 500           |
| Piece(s) / pallet             | 33000         |
| Depth of box                  | 181 mm        |
| Width of box                  | 97 mm         |
| Height of box                 | 98 mm         |
| Depth of case                 | 498 mm        |
| Width of case                 | 373 mm        |
| Height of case                | 115 mm        |
| Case volume                   | 0.0214 cbm    |
| Weight of product             | 0.0074 kg     |
| Weight of case                | 4.52 kg       |
| EAN of inner box              | 4038917275026 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Forensic swab, round, in the tube with ventilation membrane, ISO 18385, 85 mm, viscose**

EAN case

4038917275019

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**