



Product description

| | |
|---------------------|---|
| Order number | 82.1972.700 |
| Product description | Deep Well plate, Low DNA-binding, 96 well, 2.2 ml, conical wedge, without label, PCR Performance Tested, material: PP, Edged cavities, 4 piece(s)/bag |

Product characteristics

| | |
|----------------------------|--|
| Label/ Print | without label |
| Closure type | film |
| Application | collection and storage of diagnostic specimens |
| Base shape | conical wedge |
| Number of wells | 96 |
| Volume of work | 1 μ l - 2,2 ml |
| Centrifugation max (RCF) | 4500 x g |
| Centrifugation temperature | 22 |
| Centrifugation time | 5 Min |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|-----------|
| Sample volume | 2.2 ml |
| Max. volume | 2.2 ml |
| Diameter | 7.52 mm |
| Width of product | 85 mm |
| Height of product | 41.85 mm |
| Length of product | 127.76 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Purity standard | PCR Performance Tested |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 24 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 4 |
| Piece (s) / outer case | 24 |
| Piece(s) / pallet | 2304 |
| Depth of case | 283 mm |
| Width of case | 268 mm |
| Height of case | 210 mm |
| Case volume | 0.0159 cbm |
| Weight of product | 0.0868 kg |
| Weight of case | 2.43 kg |
| EAN of inner box | 4038917958042 |
| EAN case | 4038917958059 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

