



Product description

Order number	82.1972.800
Product description	Tip comb, compatible with KingFisher? Flex/Presto, Bio Sprint 96, free from DNA/DNase/RNase/PCR inhibitors, material: PP, 2 piece(s)/bag

Product characteristics

Suitable for	KingFisher? Flex/Presto, Bio Sprint 96
Autoclavable (max. °C)	121 °C

Size

Diameter	7.52 mm
Width of product	85 mm
Height of product	41.85 mm
Length of product	127.76 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard free from DNA/DNase/RNase/PCR inhibitors
Batched yes

Packaging

Minimum order qty. 20
Type of smallest subpackaging bag
Piece(s) / inner box 2
Piece (s) / outer case 20
Piece(s) / pallet 2160
Depth of case 308 mm
Width of case 223 mm
Height of case 195 mm
Case volume 0.0134 cbm
Weight of product 0.0543 kg
Weight of case 1.28 kg
EAN of inner box 4038917958080
EAN case 4038917958097

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

