



Product description

| | |
|---------------------|--|
| Order number | 86.1201 |
| Product description | Disposal bags, volume: 7 l, (LxW): 500 x 300 mm, material: PP, transparent, film thickness: 50 µm, autoclavable up to 134 °C, with print, print Biohazard, label/print: red, 50 piece(s)/bag |

Product characteristics

| | |
|------------------------|------------|
| Label/ Print | with print |
| Colour of print/label | red |
| Liquid-tight | yes |
| Autoclavable (max. °C) | 134 °C |

Size

| | |
|-------------------|--------|
| Width of product | 300 mm |
| Length of product | 500 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |

Purity & certification

| | |
|-------------------|-------------------------|
| Product category | no medical device IVD |
| Latex-free | yes |
| Safety data sheet | no |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 25000 |
| Depth of case | 348 mm |
| Width of case | 278 mm |
| Height of case | 140 mm |
| Case volume | 0.0135 cbm |
| Weight of product | 0.0141 kg |
| Weight of case | 7.3 kg |
| EAN of inner box | 4038917085755 |
| EAN case | 4038917004381 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

