



Product description

Order number	86.1202
Product description	Disposal bags, volume: 24 l, (LxW): 780 x 400 mm, material: PP, transparent, film thickness: 50 µm, autoclavable up to 134 °C, with print, print Biohazard, label/print: red, 50 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	red
Liquid-tight	yes
Autoclavable (max. °C)	134 °C

Size

Width of product	400 mm
Length of product	780 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Disposal bags, 24 l, (LxW): 780 x 400 mm, PP,
transparent, with print Biohazard

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Latex-free yes
Safety data sheet no
Batched no

Packaging

Minimum order qty. 250
Type of smallest bag
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 250
Piece(s) / pallet 9000
Depth of case 438 mm
Width of case 268 mm
Height of case 155 mm
Case volume 0.0182 cbm
Weight of product 0.0273 kg
Weight of case 7.1 kg
EAN of inner box 4038917085762
EAN case 4038917004398

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT