



## Product description

Order number	86.1203
Product description	Disposal bags, volume: 40 l, (LxW): 780 x 600 mm, material: PP, transparent, film thickness: 50 µm, autoclavable up to 134 °C, with print, print Biohazard, label/print: red, 50 piece(s)/bag

## Product characteristics

Label/ Print	with print
Colour of print/label	red
Liquid-tight	yes
Autoclavable (max. °C)	134 °C

## Size

Width of product	600 mm
Length of product	780 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Disposal bags, 40 l, (LxW): 780 x 600 mm, PP,**  
**transparent, with print Biohazard**

**Page 2**

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Latex-free yes  
Safety data sheet no  
Batched no

**Packaging**

Minimum order qty. 250  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 250  
Piece(s) / pallet 6000  
Depth of case 290 mm  
Width of case 110 mm  
Height of case 630 mm  
Case volume 0.0201 cbm  
Weight of product 0.0415 kg  
Weight of case 10.76 kg  
EAN of inner box 4038917085779  
EAN case 4038917004404

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**