



### **Product description**

Order number	82.1582.001
Product description	Micro test plate, 96 well, slip-on lid, round base, lid assembled, material: PS, transparent, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 1 piece(s)/blister

### **Product characteristics**

Label/ Print	with plastic label
Closure type	slip-on lid
Base shape	round base
Number of wells	96
Vessel type	Plate

### **Size**

Max. volume	310 $\mu$ l
Width of product	85.2 mm
Height of product	16.55 mm
Length of product	127.8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polystyrene (PS)  
Colour of product transparent

### **Purity & certification**

Product category no medical device | IVD  
Purity standard free from DNA/DNase/RNase, free from  
pyrogens/endotoxins, non-cytotoxic, sterile  
Sterilisation Gamma irradiation  
Batched yes

### **Packaging**

Minimum order qty. 50  
Type of smallest blister  
subpackaging  
Piece(s) / inner box 1  
Piece (s) / outer case 50  
Piece(s) / pallet 3200  
Depth of case 332 mm  
Width of case 297 mm  
Height of case 215 mm  
Case volume 0.0212 cbm  
Weight of product 0.0498 kg  
Weight of case 3.19 kg  
EAN of inner box 4038917131063  
EAN case 4038917347969

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

