



Product description

Order number	82.1586
Product description	Adhesive film, for ELISA and microtest plates, material: acetate, transparent, self-adhesive, 100 piece(s)/bag

Product characteristics

Suitable for	ELISA and microtest plates
--------------	----------------------------

Material & colours

Product material	acetate
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	30000
Depth of case	95 mm
Width of case	239 mm
Height of case	225 mm
Case volume	0.0051 cbm
Weight of product	0.0028 kg
Weight of case	3 kg
EAN of inner box	4038917087148
EAN case	4038917005609

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

