



### **Product description**

Order number	95.950
Product description	Felt marker, black, alcohol-resistant, application e.g.: for labelling plastic tubes and reaction tubes, 10 piece(s)/case

### **Product characteristics**

Label/ Print	with print
Colour of print/label	black
Example of use	for labelling plastic tubes and reaction tubes

### **Size**

Diameter	13 mm
Length of product	133 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	white/black

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      yes

### **Packaging**

Minimum order qty.                      10  
Type of smallest                              case  
subpackaging  
Piece(s) / inner box                      10  
Piece (s) / outer case                      500  
Piece(s) / pallet                              45000  
Depth of case                              145 mm  
Width of case                                112 mm  
Height of case                               58 mm  
Case volume                                0.0009 cbm  
Weight of product                        0.0161 kg  
Weight of case                              9.1 kg  
EAN of inner box                        4038917380355  
EAN case                                    4038917380348

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

