



### **Product description**

Order number	95.955
Product description	Felt marker, green, waterproof, application e.g.: for labelling plastic tubes and reaction tubes, 10 piece(s)/case

### **Product characteristics**

Colour of print/label	white
Example of use	for labelling plastic tubes and reaction tubes

### **Size**

Diameter	10 mm
Length of product	142 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	black

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      yes

### **Packaging**

Minimum order qty.                      10  
Type of smallest                              case  
subpackaging  
Piece(s) / inner box                      10  
Piece (s) / outer case                      1000  
Piece(s) / pallet                              48000  
Depth of case                                300 mm  
Width of case                                 291 mm  
Height of case                                392 mm  
Case volume                                 0.0342 cbm  
Weight of product                         0.009 kg  
Weight of case                                10 kg  
EAN of inner box                            4038917084642  
EAN case                                       4038917003414

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

