



### **Product description**

Order number 95.1760.100  
Product description Rubber foot, for pipette tip boxes, grey

### **Product characteristics**

Suitable for pipette tip boxes

### **Material & colours**

Product material Thermoplastic elastomers (TPE)  
Colour of product grey

### **Purity & certification**

Product category no medical device | IVD  
Batched no

### **Packaging**

Minimum order qty. 25  
Piece(s) / inner box 25  
Piece (s) / outer case 2000  
Piece(s) / pallet 648000  
Depth of case 268 mm  
Width of case 191 mm  
Height of case 100 mm  
Case volume 0.0051 cbm  
Weight of product 0.0001 kg  
Weight of case 0.45 kg  
EAN of inner box 4038917522670  
EAN case 4038917522663

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

