



Product description

Order number	70.1183.102
Product description	Pipette tip, working volume: 5 ml, transparent, filling level rings, suitable for Eppendorf, 50 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	Eppendorf
Volume of work	5 ml
Filter	no
Autoclavable (max. °C)	121 °C

Size

Max. volume	5 ml
Length of product	123 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 500
Type of smallest subpackaging box
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 16500
Depth of case 165 mm
Width of case 311 mm
Height of case 610 mm
Case volume 0.0313 cbm
Weight of product 0.0065 kg
Weight of case 3.7 kg
EAN of inner box 4038917126472
EAN case 4038917080477

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

