



### **Product description**

Order number	70.1193.150
Product description	Filter tip, working volume: 1,000 µl, black, suitable for application on SARSTEDT IVARO AP devices, 480 piece(s)/blister

### **Product characteristics**

Suitable for	application on SARSTEDT IVARO AP devices
Volume of work	1.000 µl
Filter	yes

### **Size**

Max. volume	1000 µl
Length of product	95 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	black

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE supplier
Batched	yes

### **Packaging**

Minimum order qty.	3840
Type of smallest subpackaging	blister
Piece(s) / inner box	3840
Piece (s) / outer case	3840
Piece(s) / pallet	61440
Depth of case	563 mm
Width of case	343 mm
Height of case	500 mm
Case volume	0.0966 cbm
Weight of product	0.0013 kg
Weight of case	7.9 kg
EAN of inner box	4038917434768

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Filter tip, 1,000 µl, black, 480 piece(s)/blister**

**Page 2**

EAN case

4038917434751

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**