




Product description

Order number	70.3010.110
Product description	Pipette tip, working volume: 10 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpiquette and Brand and products of identical design, 480 piece(s)/StackPack

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpiquette and Brand and products of identical design
Volume of work	10 µl
Filter	no
Coding colour	grey

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

Pipette tip, 10 µl, transparent, PCR Performance Tested, Low retention, 480 piece(s)/StackPack

Page 2

Size

Max. volume	10 µl
Length of product	31 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	1920
Type of smallest subpackaging	StackPack
Piece(s) / inner box	1920
Piece (s) / outer case	1920
Piece(s) / pallet	307200
Depth of case	290 mm
Width of case	143 mm
Height of case	225 mm
Case volume	0.0093 cbm
Weight of product	0.0004 kg
Weight of case	1.69 kg
EAN of inner box	4038917894494
EAN case	4038917894500

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT