



Product description

Order number	70.3010.210
Product description	Pipette tip, working volume: 10 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 96 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design
Volume of work	10 µl
Filter	no
Coding colour	grey
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Pipette tip, 10 µl, transparent, PCR Performance Tested, Low retention, 96 piece(s)/box

Size

Max. volume	10 µl
Length of product	31 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	1920
Type of smallest subpackaging	box
Piece(s) / inner box	480
Piece (s) / outer case	1920
Piece(s) / pallet	96000
Depth of box	330 mm
Width of box	111 mm
Height of box	148 mm
Depth of case	348 mm
Width of case	325 mm
Height of case	245 mm
Case volume	0.0277 cbm
Weight of product	0.0009 kg
Weight of case	5.16 kg
EAN of inner box	4038917894753
EAN case	4038917894760

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

