



Product description

Order number	70.3031.210
Product description	Pipette tip, working volume: 250 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 96 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design
Volume of work	250 µl
Filter	no
Coding colour	yellow
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Pipette tip, 250 µl, transparent, PCR Performance Tested, Low retention, 96 piece(s)/box

Page 2

Size

Max. volume 250 µl
Length of product 51 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard PCR Performance Tested
Batched yes

Packaging

Minimum order qty. 1920
Type of smallest subpackaging box
Piece(s) / inner box 480
Piece (s) / outer case 1920
Piece(s) / pallet 96000
Depth of box 348 mm
Width of box 129 mm
Height of box 166 mm
Depth of case 348 mm
Width of case 325 mm
Height of case 245 mm
Case volume 0.0277 cbm
Weight of product 0.0011 kg
Weight of case 5.16 kg
EAN of inner box 4038917845762
EAN case 4038917845779

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT