



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 70.3050.355  |
| Product description | Filter tip, working volume: 1,000 µl, transparent, filling level rings, Biosphere® plus, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, 96 piece(s)/SingleRefill |

### **Product characteristics**

|                |  |
|----------------|--|
| Graduation     | yes  |
| Suitable for   | SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design |
| Volume of work | 1.000 µl   |
| Filter         | yes  |
| Coding colour  | blue   |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

|                   |         |
|-------------------|---------|
| Max. volume       | 1000 µl |
| Diameter          | 9 mm    |
| Length of product | 72 mm   |

### **Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |

### **Purity & certification**

|                           |                               |
|---------------------------|-------------------------------|
| Product category          | In vitrodiagnostic device, CE |
| Certificate of conformity | yes                           |
| CE certified              | CE                            |
| Purity standard           | Biosphere® plus               |
| Sterilisation             | Gas sterilisation             |
| Pyrogen-free              | yes                           |
| Batched                   | yes                           |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1920          |
| Type of smallest subpackaging | SingleRefill  |
| Piece(s) / inner box          | 480           |
| Piece (s) / outer case        | 1920          |
| Piece(s) / pallet             | 69120         |
| Depth of box                  | 495 mm        |
| Width of box                  | 105 mm        |
| Height of box                 | 148 mm        |
| Depth of case                 | 512 mm        |
| Width of case                 | 325 mm        |
| Height of case                | 237 mm        |
| Case volume                   | 0.0394 cbm    |
| Weight of product             | 0.0015 kg     |
| Weight of case                | 4.29 kg       |
| EAN of inner box              | 4038917424608 |
| EAN case                      | 4038917424592 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

