



### **Product description**

Order number	92.1417.002
Product description	Silicon adapter, for M-Sarpette® , transparent, silicone, 1 piece(s)/bag

### **Product characteristics**

Suitable for	M-Sarpette®
--------------	-------------

### **Material & colours**

Product material	silicone
Colour of product	transparent

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	100
Piece(s) / pallet	32400
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.006 kg
Weight of case	0.85 kg
EAN of inner box	4038917129008
EAN case	4038917021906

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

